



🔍 Use a keyword, test name or number

Volatiles, Whole Blood

TEST: 007062 CPT: 80320

Test Includes Acetone; ethanol; isopropanol; methanol

Special Instructions Do not prepare venipuncture site with alcohol or remove stopper from tube.

Expected Turnaround Time 2 - 3 days

Turnaround Time Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Related Documents

- [Sample Report](#)

SPECIMEN REQUIREMENTS

Specimen Whole blood **or** serum

Volume 7 mL

Minimum Volume 0.5 mL

Container Gray-top (sodium fluoride) tube (preferred) **or** red-top tube. Submit original **unopened** tube.

Storage Instructions Room temperature

Stability Requirements	Temperature	Period
	Room temperature	14 days
	Refrigerated	14 days
	Frozen	14 days
	Freeze/thaw cycles	Stable x3

TEST DETAILS

Use Evaluate methanol and isopropanol toxicity and alcohol drug abuse

Limitations This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

Methodology Gas chromatography (GC)

© 2023 Laboratory Corporation of America® Holdings and Lexi-Comp Inc. All Rights Reserved.

[CPT Statement/Profile Statement](#)

The LOINC® codes are copyright © 1994-2023, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC) Committee. Permission is granted in perpetuity, without payment of license fees or royalties, to use, copy, or distribute the LOINC® codes for any commercial or non-commercial purpose, subject to the terms under the license agreement found at <https://loinc.org/license/>. Additional information regarding LOINC® codes can be found at LOINC.org, including the LOINC Manual, which can be downloaded at [LOINC.org/downloads/files/LOINCManual.pdf](https://loinc.org/downloads/files/LOINCManual.pdf)

