



🔍 Use a keyword, test name or number

Drug Analysis Profile, Comprehensive, Urine

TEST: 791194

CPT: 80307; 80326; 80331; 80334; 80337; 80338; 80341; 80344; 80347; 80348;
80353; 80354; 80355; 80357; 80358; 80359; 80360; 80361; 80364; 80365; 80366;
80367; 80368; 80370; 80371; 80372; 80373; 80377; 83992

Test Includes Alcohol; amphetamines; analgesics/NSAIDs; anesthetics (local); anticonvulsants; antidepressants; antihistamines; antipsychotics; barbiturates; benzodiazepines; buprenorphine; cannabinoids (THC); cocaine metabolite; creatinine; fentanyl and analogues; methadone; muscle relaxants; opiates; other hallucinogens; other opioids; oxycodone; sedatives/hypnotics; sympathomimetics; tapentadol; and miscellaneous drugs. For a complete list of profile components, see **Drug Analysis Profile, Comprehensive (791194)**, under Related Documents.

Expected 5 - 7 days

Turnaround

Time

Turnaround time is defined as the usual number of days from the date of pickup of a specimen for testing to when the result is released to the ordering provider. In some cases, additional time should be allowed for additional confirmatory or additional reflex tests. Testing schedules may vary.

Related

For more information, please view the literature below.

Documents

[Drug Analysis Profile, Comprehensive, Urine \(791194\)](#)

SPECIMEN REQUIREMENTS

Specimen

Urine (random or spot collection)

Volume

30 mL

Minimum Volume 3 mL (**Note:** This volume does **not** allow for repeat testing.)

Storage Instructions Room temperature is acceptable for shipping and/or short-term storage (three days); specimen may also be refrigerated or frozen.

Causes for Rejection Incorrect sample type; name discrepancy; no name on sample; no name on test request form

TEST DETAILS

Use Detect the presence of abused, prescribed, and over-the-counter drugs in urine

Methodology Immunoassay (IA) and liquid chromatography/tandem mass spectrometry (LC/MS-MS)

© 2023 Laboratory Corporation of America® Holdings and Lexi-Comp Inc. All Rights Reserved.

CPT Statement/Profile Statement

The LOINC® codes are copyright © 1994-2023, Regenstrief Institute, Inc. and the Logical Observation Identifiers Names and Codes (LOINC) Committee. Permission is granted in perpetuity, without payment of license fees or royalties, to use, copy, or distribute the LOINC® codes for any commercial or non-commercial purpose, subject to the terms under the license agreement found at <https://loinc.org/license/>. Additional information regarding LOINC® codes can be found at LOINC.org, including the LOINC Manual, which can be downloaded at [LOINC.org/downloads/files/LOINCManual.pdf](https://loinc.org/downloads/files/LOINCManual.pdf)

[Privacy Statement](#) [Terms and Conditions](#) [Notice of Nondiscrimination](#)

[Combatting Modern Slavery and Human Trafficking Statement](#)

© 2023 Laboratory Corporation of America® Holdings. All Rights Reserved.